



M. Stowell, T. Nemkov, Y. Peng, C. Bingham, C. Orser, K. Loh, L. Schanuel, M. Plam
AmideBio, LLC, 331 South 104th Street Louisville CO 80027, USA

Thioflavin T Assay for Fibrillation




Thioflavin T (ThT) dye fluorescence is commonly used to monitor and quantify the formation of amyloid fibrils. Interaction between ThT and β -sheets of A β fibrils is monitored by the increase in 440 nm excitation – 490 nm emission and thus increases with amyloid fibrillation. BioPure™ A β 42 reaches the maximum Δ 480 significantly faster and this rate is consequently slowed when mixed with 10% Synthetic A β 42.

MTT Assay for Cytotoxicity



Undifferentiated PC12 cells were treated with both BioPure™ A β 42 and synthetic A β 42. Cell death was monitored using the MTT assay and showed BioPure™ A β 42 to be significantly more cytotoxic than its synthetic counterpart.

Electron Micrograph of A β 42 Fibrillation



A β 42 was oligomerized in 100mM NaCl, 10mM NaH₂PO₄, 1mM HCl, pH 7.4 and incubated at room temperature for 1 hour. Samples were visualized using a Philips CM100 TEM at 80 kV. BioPure™ A β 42 (left panel) shows longer, more ordered fibrils than Synthetic A β 42 (right panel).

Amyloid- β 1-42


Amyloid- β 1-42 (A β 42) plays a central role in the pathogenesis of Alzheimer's disease. As a reagent, the purity of A β 42 is crucial for drug development.

Recombinant (BioPure™) A β 42 has reproducibly been shown to have higher *in vitro* toxicity and aggregate significantly faster than synthetically prepared A β 42 due to a higher level of purity. These traits allow for more accurate and consistent results.

BioPure™ A β 42 produced using AmideBio's proprietary technology produces material that offers the following advantages:

- >99% purity by RP-HPLC and SDS-PAGE, with mass confirmation by ESI-MS.
- Devoid of *E. coli* protein contaminants.
- Absence of (n-1) deletion products and racemized amino acids that are characteristic of chemical synthesis and highly difficult to remove.
- Batch-to-batch consistency for reproducible results.

BioPure™ A β 42 Has Higher Purity Than Synthetic A β 42



>99% BioPure™ A β 42

>95% Pure Synthetic A β 42

Most A β 42 preparations are >95% pure. AmideBio's BioPure™ A β 42 boasts a purity of >99%, which does not contain the same synthetic artifacts that have been shown to inhibit bioactivity.

Conclusions


BioPure™ A β 42 is More Biologically Active

- ✓ Faster rate of fibrillation as shown by the ThT Assay
- ✓ Higher *in vitro* cytotoxicity as shown by the MTT Assay
- ✓ Longer fibrils due to higher homogeneity

Minor Impurities Can Have Major Impacts


- ✓ Can alter activity
- ✓ Can alter formulation stability
- ✓ Cannot be easily removed from synthetic preparations

BioPure™ A β 42 Lacks (n-1) and Racemic Impurities




Comparative C18 RP-HPLC analysis of BioPure™ A β 42 to various commercially available synthetic A β 42 preparations. **A** H₂O + 0.05% TFA; **B** CH₃CN + 0.05% TFA; 60°C; 214 nm, 30-50% B over 20 minutes.

BioPure™ A β 42 has consistent batch-to-batch purity



Two separate batches of BioPure™ A β 42 produced by AmideBio and a synthetic competitor compared by ESI-MS to show both the improved purity of BioPure™ A β 42 and the batch-to-batch variation of impurities in synthetic preparations that can lead to inconsistent results.

SDS-PAGE Characterization of BioPure™ A β 42



BioPure™ A β 42 analyzed on 10-20% Tris-Glycine gel to confirm the removal of all *E. coli* protein contaminants.